Science with Passion

Azura

Click

Supplement/Beiblatt

Document No. V6711

Note: Read the corresponding technical documentation for handling and safety reasons.

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Version number: 2.1

Release date: May 02, 2023

The information in this document is subject to change without prior notice. For the latest version of the instructions, visit our website: www.knauer.net/library.

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Table of contents

1. Product information	1
2. Attaching AZURA® Click	2
3. Attaching devices on AZURA® Click	3
3.1 Interface box, pressure control, pH sensor, Organizer & distribution box	3
3.1.1 Attaching device	3
3.1.2 Taking off device	5
3.2 Air sensor	6
3.2.1 Attaching air sensor	6
3.2.2 Taking off air sensor	7
4. Repeat orders	8
5. Technical data	8

1. Product information

Fig. 1 Carrier rail AZURA® Click

The carrier rail AZURA® Click is intended for mounting the following devices:

- Interface box (AZB00XA)
- Pressure control (AZG10)
- AZURA® pH sensor (A70091-2)
- AZURA® Organizer (A70085)
- Distribution box (AZS80SA)
- Air sensor (A70092-1, A70093-1, A70083-1)

2. Attaching AZURA® Click

Prerequisite: ■ The side panel is clean, dry, and free from dust and grease.

Process

1. Loosen the knurled-head screw.
 2. Hook the left side of the carrier rail AZURA® Click into the vent holes of the device.
 3. Position the carrier rail AZURA® Click and tighten the knurled-head screw.
-

Video: On the KNAUER YouTube channel you will find a short tutorial on how to attach and take off the carrier rail: [AZURA® Click](#).

3. Attaching devices on AZURA® Click

3.1 Interface box, pressure control, pH sensor, Organizer & distribution box

3.1.1 Attaching device

Note: Interface box, pressure control, pH sensor and distribution box are supplied with a pre-installed mounting rail holder that can be used to mount the device to AZURA® Click.

Prerequisite: AZURA® Click has been mounted.

Process

1. Hook the mounting rail holder ① on the upper border of the carrier rail ②.
2. Push the device against the lower border of the carrier rail.

Fig. 2 Position of the interface box


Fig. 3 Position of the pressure control

Process

Fig. 4 Position of the AZURA® pH sensor

Fig. 5 Position of the AZURA® Organizer

Note: The mounting rail holder of the distribution box (AZS80SA) is attached to the carrier rail in the same way as for the devices in figure 4 to 7.

Video: On the KNAUER YouTube channel you will find a short tutorial on how to attach and take off the AZURA® Organizer from the carrier rail: [AZURA® Click with AZURA® Organizer](#).

3.1.2 Taking off device

Process

1. Pull the handle of the mounting rail holder down ① .
2. Take off the device by tilting it upwards ② .

Fig. 6 Taking off the device exemplified by the interface box

3.2 Air sensor

3.2.1 Attaching air sensor

- Prerequisite:**
- The metal rail has been removed (Fig. 7).
 - AZURA® Click has been mounted.

Process

1. Pull back the metal rail ① against the spring.
2. Take off the metal rail of the air sensor.

Fig. 7 Removing the metal rail of the air sensor

Process

Variant 1:

3. Hook the lower side of the air sensor onto the lower border of the carrier rail ②.
4. Pull the air sensor slightly upwards.
5. Push the upper side of the device against the upper border of the carrier rail until the air sensor clicks into place.

Variant 2:

3. Hook the upper side of the air sensor onto the upper border of the carrier rail.
4. Pull the air sensor slightly downwards.
5. Push the lower side of the device against the lower border of the carrier rail until the air sensor clicks into place.

Fig. 8 Position of air sensor

Video: On the KNAUER YouTube channel you will find a short tutorial on how to attach and take off the air sensor from the carrier rail: [AZURA® Click with air sensor](#).

3.2.2 Taking off air sensor

Process

1. Pull the air sensor over the top of the upper carrier rail.
2. Take off the air sensor.

4. Repeat orders

Name	Article no.
AZURA® Click	A70089
Interface box	AZB00XA
Air sensor for 1/16" tubing	A70092-1
Air sensor for 1/8" tubing	A70093-1
Air sensor for 1/4" tubing	A70083-1
AZURA® pH sensor	A70091-2
Pressure control	AZG10
Distribution box	AZS80SA

5. Technical data

Dimensions (B × H × T)	330 × 35 × 7.5 mm
Material	Stainless steel

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All rights reserved.
AZURA® is a registered trademark of KNAUER
Wissenschaftliche Geräte GmbH